

PORTABLE MINE ALNU X-RAY




40~80kV
3mA
0.4mm
FOCAL SPOT

HIGH- RESOLUTION IMAGES

IMAGE QUALITY IS IMPROVED
BY REDUCTION OF RIPPLES AND PHANTOM PHENOMENON
THROUGH HIGH FREQUENCY INVERTER.

m-SENSOR
17 x 17 DETECTOR



SPECIFICATION

X-ray Generator

TYPE	High Frequency Inverter
Power Output	120W ~ 240W
Tube Voltage	40 ~ 80 kV(Increment : 1kV)
Tube Current	3 mA (Fixed)
Exposure Time Range	0.1 sec ~ 1.3 sec (Increment : 0.1 sec)
Focal Spot Size	0.4mm

Power Adapter

Input	AC 100-240V, 50/60Hz, 1.0A
Output	DC 12V, 2.1A(MAX)

Battery

TYPE	Lithium-ion Battery
Battery Capacity	2500mAh
Supply Voltage	DC 25.2V